




ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

BI-LUMEN IRRIGATION EMBOLECTOMY CATHETER (DISTAL / PROXIMAL)



ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET / SOFT TIP / RADIOPAQUE

INTENDED FOR INTRA-OPERATIVE REMOVAL OF ARTERIAL EMBOLI FROM ARTERIES

AVAILABLE SIZES:

MODEL: DISTAL AND PROXIMAL

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
4F	RED	4F IRR 40 - DISTAL / PROXIMAL 4F IRR 80 - DISTAL / PROXIMAL	0.50 CC	9 MM
5F	WHITE	5F IRR 40 - DISTAL / PROXIMAL 5F IRR 80 - DISTAL / PROXIMAL	0.75 CC	11 MM
6F	BLUE	6F IRR 40 - DISTAL / PROXIMAL 6F IRR 80 - DISTAL / PROXIMAL	1.25 CC	13 MM
7F	YELLOW	7F IRR 40 - DISTAL / PROXIMAL 7F IRR 80 - DISTAL / PROXIMAL	2.00 CC	14 MM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: lmi@earthlink.net

www.lucasmedicalinc.com

ISO 13485:2012 & ISO 9001-2008 REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

MBM
C/CERIFALTE 117, 28230 LAS ROZAS
MADRID, SPAIN
TEL: +34916308273



0
4
7
3