



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

## BILIARY CATHETER



ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET  
SOFT TIP / RADIOPAQUE

*INTENDED FOR REMOVAL OF GALL STONES FROM THE BILIARY DUCT*

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
3F	GREEN	3F BIL 23	0.10 CC	5 MM
3F	GREEN	3F BIL 40	0.10 CC	5 MM
4F	RED	4F BIL 23	0.50 CC	9 MM
4F	RED	4F BIL 40	0.50 CC	9 MM
5F	WHITE	5F BIL 23	0.75 CC	11 MM
5F	WHITE	5F BIL 40	0.75 CC	11 MM
6F	BLUE	6F BIL 23	1.25 CC	13 MM
6F	BLUE	6F BIL 40	1.25 CC	13 MM
7F	YELLOW	7F BIL 23	2.00 CC	14 MM
7F	YELLOW	7F BIL 40	2.00 CC	14 MM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: [lmi@earthlink.net](mailto:lmi@earthlink.net)

[www.lucasmedicalinc.com](http://www.lucasmedicalinc.com)

ISO 13485:2012 & ISO 9001-2008 REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

MBM  
C/CERIFALTE 117, 28230 LAS ROZAS  
MADRID, SPAIN  
TEL: +34916308273



0  
4  
7  
3