




ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

CAROTID ENDARTERECTOMY SHUNT CATHETER

INTENDED FOR PROVIDING TEMPORARY BLOOD FLOW TO CEREBRAL AREA DURING A CAROTID ENDARTERECTOMY PROCEDURE



- ❖ ALL SILICONE
- ❖ NON-FRAGMENTING BALLOON
- ❖ SOFT TIP
- ❖ RADIOPAQUE

AVAILABLE SIZES:

MODEL: C E S B

FRENCH SIZE	PRODUCT DESIGNATION	RECOMMENDED LIQUID VOLUME (CC)		INFLATED DIAMETER (MM)		LENGTH (CM)
		INTERNAL BALLOON	COMMON BALLOON	INTERNAL BALLOON	COMMON BALLOON	
8F	8 CESB 15	3/4	1 1/4	11	12	15
8F	8 CESB 30	3/4	1 1/4	11	12	30
9F	9 CESB 15	3/4	1 1/4	12	13	15
9F	9 CESB 30	3/4	1 1/4	12	13	30
10F	10 CESB 15	3/4	1 1/4	12	13.5	15
10F	10 CESB 30	3/4	1 1/4	12	13.5	30
11F	11 CESB 15	3/4	1 1/4	12.5	13.5	15
11F	11 CESB 30	3/4	1 1/4	12.5	13.5	30
12F	12 CESB 15	3/4	1 1/4	12.5	14	15
12F	12 CESB 30	3/4	1 1/4	12.5	14	30

1751 S. DOUGLASS RD, ANAHEIM, CA 92806
E-MAIL: lmi@earthlink.net

PHONE: (714) 938-0233 FAX: (714) 938-0130
www.lucasmedicalinc.com

ISO 13485:2012 & ISO 9001-2008 REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

MBM
C/CERIFALTE 117, 28230 LAS ROZAS
MADRID, SPAIN
TEL: +34916308273

