



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

ARTERIAL EMBOLECTOMY CATHETER



**ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET
SOFT TIP / RADIOPAQUE**

INTENDED FOR INTRA-OPERATIVE REMOVAL OF ARTERIAL EMBOLI FROM ARTERIES

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
2F	PURPLE	2F EMB 40	0.05 CC	4 MM
2F	PURPLE	2F EMB 60	0.05 CC	4 MM
3F	GREEN	3F EMB 40	0.10 CC	5 MM
3F	GREEN	3F EMB 80	0.10 CC	5 MM
4F	RED	4F EMB 40	0.50 CC	9 MM
4F	RED	4F EMB 80	0.50 CC	9 MM
4F	RED	4F EMB 100	0.50 CC	9 MM
5F	WHITE	5F EMB 40	0.75 CC	11 MM
5F	WHITE	5F EMB 80	0.75 CC	11 MM
6F	BLUE	6F EMB 40	1.25 CC	13 MM
6F	BLUE	6F EMB 80	1.25 CC	13 MM
7F	YELLOW	7F EMB 40	2.00 CC	14 MM
7F	YELLOW	7F EMB 80	2.00 CC	14 MM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: lmi@earthlink.net

www.lucasmedicalinc.com

ISO 13485:2012 & ISO 9001-2008 REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

MBM
C/CERIFALTE 117, 28230 LAS ROZAS
MADRID, SPAIN
TEL: +34916308273

