



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

OCCLUSION CATHETER



ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET
SOFT TIP / RADIOPAQUE

INTENDED FOR SHORT TERM OCCLUSION OF BLOOD FLOW

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	MAXIMUM LIQUID INFLATION VOLUME (MLV)	INFLATED BALLOON DIAMETER @ MLV
3F	GREEN	3F Occ 23	0.10 CC	5 MM
3F	GREEN	3F Occ 40	0.10 CC	5 MM
4F	RED	4F Occ 23	0.50 CC	9 MM
4F	RED	4F Occ 40	0.50 CC	9 MM
5F	WHITE	5F Occ 23	0.75 CC	11 MM
5F	WHITE	5F Occ 40	0.75 CC	11 MM
6F	BLUE	6F Occ 23	1.25 CC	13 MM
6F	BLUE	6F Occ 40	1.25 CC	13 MM
7F	YELLOW	7F Occ 23	2.00 CC	14 MM
7F	YELLOW	7F Occ 40	2.00 CC	14 MM

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: lmi@earthlink.net

www.lucasmedicalinc.com

ISO 13485:2012 & ISO 9001-2008 REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

MBM
C/CERIFALTE 117, 28230 LAS ROZAS
MADRID, SPAIN
TEL: +34916308273



0
4
7
3