



ADVANCED TECHNOLOGY

ADDRESSING TODAY'S NEEDS

LATEX FREE

IRRIGATION CATHETER



**ALL SILICONE / NON-FRAGMENTING BALLOON / INTEGRATED FLEXIBLE STYLET
SOFT TIP / RADIOPAQUE**

*INTENDED FOR HEPARINIZING THE ARTERIAL TREE AND VEINS FOLLOWING EMBOLECTOMY AND
THROMBECTOMY PROCEDURES AS WELL AS THE BILIARY TREE AFTER CHOLANGIOGRAPHY*

AVAILABLE SIZES:

FRENCH SIZE	COLOR CODE	PRODUCT DESIGNATION	LENGTH (CM)
4F	RED	4F IGR 80	80
5F	WHITE	5F IGR 80	80
6F	BLUE	6F IGR 80	80
7F	YELLOW	7F IGR 80	80

1751 S. DOUGLASS RD, ANAHEIM, CA 92806

PHONE: (714) 938-0233 FAX: (714) 938-0130

E-MAIL: lmi@earthlink.net

www.lucasmedicalinc.com

ISO 13485:2012 & ISO 9001-2008 REGISTERED

U.S. PATENT NUMBERS: 5,762,996 & 5,795,332

MBM
C/CERIFALTE 117, 28230 LAS ROZAS
MADRID, SPAIN
TEL: +34916308273

